

# RASAL1 Antibody, HRP conjugated



PACO56459

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable inhibitory regulator of the Ras-cyclic AMP pathway. Plays a role in dendrite formation by melanocytes.

**Gene ID:**

RASAL1

**Uniprot**

O95294

**Synonyms:**

RasGAP-activating-like protein 1 (RAS protein activator like 1) (Ras GTPase-activating-like protein), RASAL1, RASAL

**Immunogen:**

Recombinant Human RasGAP-activating-like protein 1 protein (698-804AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A