

RASAL1 Antibody, FITC conjugated



PACO56460

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable inhibitory regulator of the Ras-cyclic AMP pathway. Plays a role in dendrite formation by melanocytes.

Gene ID:

RASAL1

Uniprot

O95294

Synonyms:

RasGAP-activating-like protein 1 (RAS protein activator like 1) (Ras GTPase-activating-like protein), RASAL1, RASAL

Immunogen:

Recombinant Human RasGAP-activating-like protein 1 protein (698-804AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A