

# ASMT Antibody, FITC conjugated



PACO56730

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Isoform 1 catalyzes the transfer of a methyl group onto N-acetylserotonin, producing melatonin (N-acetyl-5-methoxytryptamine). Isoform 2 and isoform 3 lack enzyme activity.

**Gene ID:**

ASMT

**Uniprot**

P46597

**Synonyms:**

Acetylserotonin O-methyltransferase (EC 2.1.1.4) (Hydroxyindole O-methyltransferase) (HIOMT), ASMT

**Immunogen:**

Recombinant Human Acetylserotonin O-methyltransferase protein (221-316AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A