

lacZ Antibody



PACO56804

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Hydrolysis of terminal non-reducing beta-D-galactose residues in beta-D-galactosides.

Gene ID:

lacZ

Uniprot

P00722

Synonyms:

Beta-galactosidase (Beta-gal) (EC 3.2.1.23) (Lactase), lacZ


Immunogen:

Recombinant Escherichia coli Beta-galactosidase protein (2-1024AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Recombinant protein. All lanes: lacZ antibody at 7.5 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 117 kDa. Observed band size: 117 kDa.