

# GPR1 Antibody, HRP conjugated



PACO56809

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Receptor for the inflammation-associated leukocyte chemoattractant chemerin/RARRES2 suggesting a role for this receptor in the regulation of inflammation. Can act as a coreceptor for HIV-1.

**Gene ID:**

GPR1

**Uniprot**

P46091

**Synonyms:**

G-protein coupled receptor 1, GPR1

**Immunogen:**

Recombinant Human G-protein coupled receptor 1 protein (308-355AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A