

HOXA4 Antibody, HRP conjugated



PACO56837

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Binds to sites in the 5'-flanking sequence of its coding region with various affinities. The consensus sequences of the high and low affinity binding sites are 5'-TAATGA[CG]-3' and 5'-CTAATTTT-3'.

Gene ID:

HOXA4

Uniprot

Q00056

Synonyms:

Homeobox protein Hox-A4 (Homeobox protein Hox-1.4) (Homeobox protein Hox-1D), HOXA4, HOX1D

Immunogen:

Recombinant Human Homeobox protein Hox-A4 protein (23-173AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A