

HOXD3 Antibody



PACO56844

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Gene ID:

HOXD3

Uniprot

P31249

Synonyms:

Homeobox protein Hox-D3 (Homeobox protein Hox-4A), HOXD3, HOX1D HOX4A


Immunogen:

Recombinant Human Homeobox protein Hox-D3 protein (263-346AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Jurkat whole cell lysate, K562 whole cell lysate, SH-SY5Y whole cell lysate, Rat brain tissue. All lanes: HOXD3 antibody at 8.5 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 46 kDa. Observed band size: 46 kDa.