

HOXD8 Antibody



PACO56848

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Gene ID:

HOXD8

Uniprot

P13378

Synonyms:

Homeobox protein Hox-D8 (Homeobox protein Hox-4E) (Homeobox protein Hox-5.4), HOXD8, HOX4E


Immunogen:

Recombinant Human Homeobox protein Hox-D8 protein (59-108AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: PC-3 whole cell lysate, Mouse brain tissue. All lanes: HOXD8 antibody at 4.5 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 32, 9 kDa. Observed band size: 32 kDa.