

LHX2 Antibody, HRP conjugated



PACO56873

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acts as a transcriptional activator. Stimulates the promoter of the alpha-glycoprotein gene. Transcriptional regulatory protein involved in the control of cell differentiation in developing lymphoid and neural cell types.

Gene ID:

LHX2

Uniprot

P50458

Synonyms:

LIM/homeobox protein Lhx2 (Homeobox protein LH-2) (LIM homeobox protein 2), LHX2, LH2

Immunogen:

Recombinant Human LIM/homeobox protein Lhx2 protein (197-377AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A