

# LRP2 Antibody, HRP conjugated



PACO56881

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Acts together with cubilin to mediate HDL endocytosis. May participate in regulation of parathyroid-hormone and para-thyroid-hormone-related protein release.

**Gene ID:**

LRP2

**Uniprot**

P98164

**Synonyms:**

Low-density lipoprotein receptor-related protein 2 (LRP-2) (Glycoprotein 330) (gp330) (Megalin), LRP2

**Immunogen:**

Recombinant Human Low-density lipoprotein receptor-related protein 2 protein (1211-1369AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A