

PHOX2A Antibody, HRP conjugated



PACO56925

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in regulating the specificity of expression of the catecholamine biosynthetic genes. Acts as a transcription activator/factor. Could maintain the noradrenergic phenotype.

Gene ID:

PHOX2A

Uniprot

O14813

Synonyms:

Paired mesoderm homeobox protein 2A (ARX1 homeodomain protein) (Aristaless homeobox protein homolog) (Paired-like homeobox 2A), PHOX2A, ARX PMX2A

Immunogen:

Recombinant Human Paired mesoderm homeobox protein 2A protein (150-264AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A