

# SERPINB4 Antibody, Biotin conjugated



PACO56975

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May act as a protease inhibitor to modulate the host immune response against tumor cells.

**Gene ID:**

SERPINB4

**Uniprot**

P48594

**Synonyms:**

Serpin B4 (Leupin) (Peptidase inhibitor 11) (PI-11) (Squamous cell carcinoma antigen 2) (SCCA-2), SERPINB4, PI11 SCCA2

**Immunogen:**

Recombinant Human Serpin B4 protein (128-208AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A