

# KIAA0319 Antibody, HRP conjugated



PACO57045

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in neuronal migration during development of the cerebral neocortex. May function in a cell autonomous and a non-cell autonomous manner and play a role in appropriate adhesion between migrating neurons and radial glial fibers. May also regulate growth and differentiation of dendrites.

**Gene ID:**

KIAA0319

**Uniprot**

Q5VV43

**Synonyms:**

Dyslexia-associated protein KIAA0319, KIAA0319

**Immunogen:**

Recombinant Human Dyslexia-associated protein KIAA0319 protein (153-268AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A