

KIAA0319 Antibody, Biotin conjugated



PACO57047

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in neuronal migration during development of the cerebral neocortex. May function in a cell autonomous and a non-cell autonomous manner and play a role in appropriate adhesion between migrating neurons and radial glial fibers. May also regulate growth and differentiation of dendrites.

Gene ID:

KIAA0319

Uniprot

Q5VV43

Synonyms:

Dyslexia-associated protein KIAA0319, KIAA0319

Immunogen:

Recombinant Human Dyslexia-associated protein KIAA0319 protein (153-268AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A