

PON3 Antibody, HRP conjugated



PACO57105

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Has low activity towards the organophosphate paraxon and aromatic carboxylic acid, esters. Rapidly hydrolyzes lactones such as statin prodrugs (e. g. lovastatin). Hydrolyzes aromatic lactones and 5- or 6-member ring lactones with aliphatic substituents but not simple lactones or those with polar substituents.

Gene ID:

PON3

Uniprot

Q15166

Synonyms:

Serum paraoxonase/lactonase 3 (EC 3.1.1.2) (EC 3.1.1.81) (EC 3.1.8.1), PON3

Immunogen:

Recombinant Human Serum paraoxonase/lactonase 3 protein (269-335AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A