

# PAFAH1B3 Antibody, HRP conjugated



PACO57133

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Inactivates paf by removing the acetyl group at the sn-2 position. This is a catalytic subunit. Plays an important role during the development of brain.

**Gene ID:**

PAFAH1B3

**Uniprot**

Q15102

**Synonyms:**

Platelet-activating factor acetylhydrolase IB subunit gamma (EC 3.1.1.47) (PAF acetylhydrolase 29 kDa subunit) (PAF-AH 29 kDa subunit) (PAF-AH subunit gamma) (PAFAH subunit gamma), PAFAH1B3, PAFAHG

**Immunogen:**

Recombinant Human Platelet-activating factor acetylhydrolase IB subunit &gamma; protein (103-187AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A