

# SIM2 Antibody, HRP conjugated



PACO57145

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcription factor that may be a master gene of CNS development in cooperation with Arnt. It may have pleiotropic effects in the tissues expressed during development.

**Gene ID:**

SIM2

**Uniprot**

Q14190

**Synonyms:**

Single-minded homolog 2 (Class E basic helix-loop-helix protein 15) (bHLHe15), SIM2, BHLHE15

**Immunogen:**

Recombinant Human Single-minded homolog 2 protein (347-433AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A