

SIM2 Antibody, Biotin conjugated



PACO57147

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcription factor that may be a master gene of CNS development in cooperation with Arnt. It may have pleiotropic effects in the tissues expressed during development.

Gene ID:

SIM2

Uniprot

Q14190

Synonyms:

Single-minded homolog 2 (Class E basic helix-loop-helix protein 15) (bHLHe15), SIM2, BHLHE15

Immunogen:

Recombinant Human Single-minded homolog 2 protein (347-433AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A