

# SHROOM2 Antibody, FITC conjugated



PACO57154

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in endothelial cell morphology changes during cell spreading. In the retinal pigment epithelium, may regulate the biogenesis of melanosomes and promote their association with the apical cell surface by inducing gamma-tubulin redistribution.

**Gene ID:**

SHROOM2

**Uniprot**

Q13796

**Synonyms:**

Protein Shroom2 (Apical-like protein) (Protein APXL), SHROOM2, APXL

**Immunogen:**

Recombinant Human Protein Shroom2 protein (213-405AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A