

KCNT1 Antibody, Biotin conjugated



PACO57179

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Outwardly rectifying potassium channel subunit that may coassemble with other Slo-type channel subunits. Activated by high intracellular sodium or chloride levels. Activated upon stimulation of G-protein coupled receptors, such as CHRM1 and GRIA1. May be regulated by calcium in the absence of sodium ions (in vitro).

Gene ID:

KCNT1

Uniprot

Q5JUK3

Synonyms:

Potassium channel subfamily T member 1 (KCa4.1), KCNT1, KIAA1422

Immunogen:

Recombinant Human Potassium channel subfamily T member 1 protein (1026-1230AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A