

UMPS Antibody, HRP conjugated



PACO57205

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, cytosol, nucleus, orotate phosphoribosyltransferase activity, orotidine-5'-phosphate decarboxylase activity, pyrimidine nucleoside biosynthetic process, UMP biosynthetic process.

Gene ID:

UMPS

Uniprot

P11172

Synonyms:

Uridine 5'-monophosphate synthase (UMP synthase) [Includes: Orotate phosphoribosyltransferase (OPRT) (OPRTase) (EC 2.4.2.10); Orotidine 5'-phosphate decarboxylase (ODC) (EC 4.1.1.23) (OMPdecase)], UMPS

Immunogen:

Recombinant Human Uridine 5'-monophosphate synthase protein (314-444AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A