

CD300LG Antibody, FITC conjugated



PACO57242

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor which may mediate L-selectin-dependent lymphocyte rollings. Binds SELL in a calcium dependent manner. Binds lymphocyte.

Gene ID:

CD300LG

Uniprot

Q6UXG3

Synonyms:

CMRF35-like molecule 9 (CLM-9) (CD300 antigen-like family member G) (Triggering receptor expressed on myeloid cells 4) (TREM-4) (CD antigen CD300g), CD300LG, CLM9 TREM4

Immunogen:

Recombinant Human CMRF35-like molecule 9 protein (128-241AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A