

# CD99L2 Antibody, Biotin conjugated



PACO57351

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Plays a role in a late step of leukocyte extravasation helping cells to overcome the endothelial basement membrane. Acts at the same site as, but independently of, PECAM1. Homophilic adhesion molecule, but these interactions may not be required for cell aggregation.

**Gene ID:**

CD99L2

**Uniprot**

Q8TCZ2

**Synonyms:**

CD99 antigen-like protein 2 (MIC2-like protein 1) (CD antigen CD99), CD99L2, MIC2L1

**Immunogen:**

Recombinant Human CD99 antigen-like protein 2 protein (26-145AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A