

GBP4 Antibody, Biotin conjugated



PACO57399

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds GTP, GDP and GMP. Hydrolyzes GTP very efficiently; GDP rather than GMP is the major reaction product. Plays a role in erythroid differentiation.

Gene ID:

GBP4

Uniprot

Q96PP9

Synonyms:

Guanylate-binding protein 4 (EC 3.6.5) (GTP-binding protein 4) (GBP-4) (Guanine nucleotide-binding protein 4), GBP4

Immunogen:

Recombinant Human Guanylate-binding protein 4 protein (494-603AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A