

# ALS2CR12 Antibody, HRP conjugated



PACO57405

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

ALS2CR12

**Uniprot**

Q96Q35

**Synonyms:**

Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 12 protein, ALS2CR12

**Immunogen:**

Recombinant Human Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 12 protein (26-76AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A