

ALS2CR12 Antibody, Biotin conjugated



PACO57407

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

ALS2CR12

Uniprot

Q96Q35

Synonyms:

Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 12 protein, ALS2CR12

Immunogen:

Recombinant Human Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 12 protein (26-76AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A