

LEMD2 Antibody, HRP conjugated



PACO57409

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in nuclear structure organization. Required for maintaining the integrity of the nuclear envelope.

Gene ID:

LEMD2

Uniprot

Q8NC56

Synonyms:

LEM domain-containing protein 2 (hLEM2), LEMD2

Immunogen:

Recombinant Human LEM domain-containing protein 2 protein (17-204AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A