

TNFRSF19 Antibody, FITC conjugated



PACO57506

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Can mediate activation of JNK and NF-kappa-B. May promote caspase-independent cell death.

Gene ID:

TNFRSF19

Uniprot

Q9NS68

Synonyms:

Tumor necrosis factor receptor superfamily member 19 (TRADE) (Toxicity and JNK inducer), TNFRSF19, TAJ TROY

Immunogen:

Recombinant Human Tumor necrosis factor receptor superfamily member 19 protein (276-415AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A