

ADAMTS13 Antibody, Biotin conjugated



PACO57655

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cleaves the vWF multimers in plasma into smaller forms thereby controlling vWF-mediated platelet thrombus formation.

Gene ID:

ADAMTS13

Uniprot

Q76LX8

Synonyms:

A disintegrin and metalloproteinase with thrombospondin motifs 13 (ADAM-TS 13) (ADAM-TS13) (ADAMTS-13) (EC 3.4.24.87) (von Willebrand factor-cleaving protease) (vWF-CP) (vWF-cleaving protease), ADAMTS13, C9orf8

Immunogen:

Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 13 protein (738-905AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A