

BCKDHA Antibody, HRP conjugated



PACO57673

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The branched-chain alpha-keto dehydrogenase complex catalyzes the overall conversion of alpha-keto acid, to acyl-CoA and CO₂. It contains multiple copies of three enzymatic components: branched-chain alpha-keto acid, decarboxylase (E1), lipoamide acyltransferase (E2) and lipoamide dehydrogenase (E3).

Gene ID:

BCKDHA

Uniprot

P12694

Synonyms:

2-oxoisovalerate dehydrogenase subunit alpha, mitochondrial (EC 1.2.4.4) (Branched-chain alpha-keto acid, dehydrogenase E1 component alpha chain) (BCKDE1A) (BCKDH E1-alpha), BCKDHA

Immunogen:

Recombinant Human 2-oxoisovalerate dehydrogenase subunit alpha, mitochondrial protein (268-407AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A