

CSN2 Antibody, FITC conjugated



PACO57726

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Important role in determination of the surface properties of the casein micelles.

Gene ID:

CSN2

Uniprot

P05814

Synonyms:

Beta-casein, CSN2, CASB

Immunogen:

Recombinant Human Beta-casein protein (16-226AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A