

ECM1 Antibody, Biotin conjugated



PACO57743

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in endochondral bone formation as negative regulator of bone mineralization. Stimulates the proliferation of endothelial cells and promotes angiogenesis. Inhibits MMP9 proteolytic activity.

Gene ID:

ECM1

Uniprot

Q16610

Synonyms:

Extracellular matrix protein 1 (Secretory component p85), ECM1

Immunogen:

Recombinant Human Extracellular matrix protein 1 protein (386-540AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A