

# GBX2 Antibody, FITC conjugated



PACO57774

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May act as a transcription factor for cell pluripotency and differentiation in the embryo.

**Gene ID:**

GBX2

**Uniprot**

P52951

**Synonyms:**

Homeobox protein GBX-2 (Gastrulation and brain-specific homeobox protein 2), GBX2

**Immunogen:**

Recombinant Human Homeobox protein GBX-2 protein (65-248AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A