

# GSS Antibody, Biotin conjugated



PACO57795

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

cytosol, extracellular exosome, ATP binding, glutathione binding, glutathione synthase activity, identical protein binding, magnesium ion binding, protein homodimerization activity, cellular amino acid, metabolic process, glutathione biosynthetic process.

**Gene ID:**

GSS

**Uniprot**

P48637

**Synonyms:**

Glutathione synthetase (GSH synthetase) (GSH-S) (EC 6.3.2.3) (Glutathione synthase), GSS

**Immunogen:**

Recombinant Human Glutathione synthetase protein (290-418AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A