

MVD Antibody, Biotin conjugated



PACO57871

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Performs the first committed step in the biosynthesis of isoprenes.

Gene ID:

MVD

Uniprot

P53602

Synonyms:

Diphosphomevalonate decarboxylase (EC 4.1.1.33) (Mevalonate (diphospho)decarboxylase) (MDDase) (Mevalonate pyrophosphate decarboxylase), MVD, MPD

Immunogen:

Recombinant Human Diphosphomevalonate decarboxylase protein (69-221AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A