

NKX2-8 Antibody, HRP conjugated



PACO57893

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

double-stranded DNA binding, RNA polymerase II core promoter proximal region sequence-specific DNA binding, sequence-specific DNA binding, transcription factor activity, sequence-specific DNA binding, transcriptional activator activity, RNA polymerase II core promoter proximal region sequence-specific binding, liver development, positive regulation of transcription from RNA polymerase II promoter, transcription from RNA polymerase II promoter, transcription, DNA-templated.

Gene ID:

NKX2-8

Uniprot

O15522

Synonyms:

Homeobox protein Nkx-2.8 (Homeobox protein NK-2 homolog H), NKX2-8, NKX-2.8 NKX2G NKX2H

Immunogen:

Recombinant Human Homeobox protein Nkx-2.8 protein (1-160AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A