

PCCB Antibody, FITC conjugated



PACO57914

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, mitochondrial matrix, mitochondrion, CoA carboxylase activity, propionyl-CoA carboxylase activity, biotin metabolic process, fatty acid, catabolic process, short-chain fatty acid, catabolic process.

Gene ID:

PCCB

Uniprot

P05166

Synonyms:

Propionyl-CoA carboxylase beta chain, mitochondrial (PCCase subunit beta) (EC 6.4.1.3) (Propanoyl-CoA: carbon dioxide ligase subunit beta), PCCB

Immunogen:

Recombinant Human Propionyl-CoA carboxylase β chain, mitochondrial protein (200-314AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A