

PCSK7 Antibody, HRP conjugated



PACO57917

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serine endoprotease that processes various proproteins by cleavage at paired basic amino acid, , recognizing the RXXX[KR]R consensus motif. Likely functions in the constitutive secretory pathway.

Gene ID:

PCSK7

Uniprot

Q16549

Synonyms:

Proprotein convertase subtilisin/kexin type 7 (EC 3.4.21) (Lymphoma proprotein convertase) (Prohormone convertase 7) (Proprotein convertase 7) (PC7) (Proprotein convertase 8) (PC8) (hPC8) (Subtilisin/kexin-like protease PC7), PCSK7, LPC PC7 PC8 SPC7

Immunogen:

Recombinant Human Proprotein convertase subtilisin/kexin type 7 protein (189-362AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A