

SLC26A4 Antibody, HRP conjugated



PACO57973

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Sodium-independent transporter of chloride and iodide.

Gene ID:

SLC26A4

Uniprot

O43511

Synonyms:

Pendrin (Sodium-independent chloride/iodide transporter) (Solute carrier family 26 member 4), SLC26A4, PDS

Immunogen:

Recombinant Human Pendrin protein (592-655AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A