

SLC7A1 Antibody, HRP conjugated



PACO57981

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

High-affinity, low capacity permease involved in the transport of the cationic amino acid, (arginine, lysine and ornithine) in non-hepatic tissues.

Gene ID:

SLC7A1

Uniprot

P30825

Synonyms:

High affinity cationic amino acid, transporter 1 (CAT-1) (CAT1) (Ecotropic retroviral leukemia receptor homolog) (Ecotropic retrovirus receptor homolog) (Solute carrier family 7 member 1) (System Y+ basic amino acid, transporter), SLC7A1, ATRC1 ERR REC1L

Immunogen:

Recombinant Human High affinity cationic amino acid, transporter 1 protein (430-492AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A