

TMSB4X Antibody



PACO58004

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Plays an important role in the organization of the cytoskeleton. Binds to and sequesters actin monomers (G actin) and therefore inhibits actin polymerization.

Gene ID:

TMSB4X

Uniprot

P62328

Synonyms:

Thymosin beta-4 (T beta-4) (Fx) [Cleaved into: Hematopoietic system regulatory peptide (Seraspenide)], TMSB4X, TB4X THYB4 TMSB4


Immunogen:

Recombinant Human Thymosin beta-4 protein (6-44AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lymph node tissue using PACO58004 at dilution of 1:100.