

UGT1A6 Antibody, FITC conjugated



PACO58026

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

UDPGT is of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds. This isoform has specificity for phenols. Isoform 3 lacks transferase activity but acts as a negative regulator of isoform 1.

Gene ID:

UGT1A6

Uniprot

P19224

Synonyms:

UDP-glucuronosyltransferase 1-6 (UDPGT 1-6) (UGT1*6) (UGT1-06) (UGT1.6) (EC 2.4.1.17) (Phenol-metabolizing UDP-glucuronosyltransferase) (UDP-glucuronosyltransferase 1-F) (UGT-1F) (UGT1F) (UDP-glucuronosyltransferase 1A6), UGT1A6, GNT1 UGT1

Immunogen:

Recombinant Human UDP-glucuronosyltransferase 1-6 protein (278-410AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A