

Beta-mammal toxin Css4 Antibody



PACO58092

Product Information

Size:

50ug

Reactivity:

Centruroides suffusus suffusus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

beta toxins bind voltage-independently at site-4 of sodium channels (Nav) and shift the voltage of activation toward more negative potentials thereby affecting sodium channel activation and promoting spontaneous and repetitive firing. This toxin is active only on mammals.

Gene ID:**Uniprot**

P60266

Synonyms:

Beta-mammal toxin Css4 (Css IV) (CcssIV)

Immunogen:

Recombinant Centruroides suffusus suffusus Beta-mammal toxin Css4 protein (1-66AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A