

AANAT Antibody, HRP conjugated



PACO58173

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Controls the night/day rhythm of melatonin production in the pineal gland. Catalyzes the N-acetylation of serotonin into N-acetylserotonin, the penultimate step in the synthesis of melatonin.

Gene ID:

AANAT

Uniprot

Q16613

Synonyms:

Serotonin N-acetyltransferase (Serotonin acetylase) (EC 2.3.1.87) (Aralkylamine N-acetyltransferase) (AA-NAT), AANAT, SNAT

Immunogen:

Recombinant Human Serotonin N-acetyltransferase protein (1-207AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A