

# PAG2 Antibody, FITC conjugated



PACO58182

---

## Product Information

**Size:**

50ug

**Reactivity:**

Bovine

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

PAG2 or a processed derivative of this molecule might represent a factor that binds the LH receptor.

**Gene ID:**

PAG2

**Uniprot**

Q28057

**Synonyms:**

Pregnancy-associated glycoprotein 2 (PAG 2) (EC 3.4.23), PAG2

**Immunogen:**

Recombinant Bovine Pregnancy-associated glycoprotein 2 protein (22-376AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A