

NLRC5 Antibody, HRP conjugated



PACO58233

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable regulator of the NF-kappa-B and type I interferon signaling pathways. May also regulate the type II interferon signaling pathway. Plays a role in homeostatic control of innate immunity and in antiviral defense mechanisms.

Gene ID:

NLRC5

Uniprot

Q86WI3

Synonyms:

Protein NLRC5 (Caterpillar protein 16.1) (CLR16.1) (Nucleotide-binding oligomerization domain protein 27) (Nucleotide-binding oligomerization domain protein 4), NLRC5, NOD27 NOD4

Immunogen:

Recombinant Human Protein NLRC5 protein (95-208AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A