

DIO2 Antibody, Biotin conjugated



PACO58295

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Responsible for the deiodination of T4 (3,5,3',5'-tetraiodothyronine) into T3 (3,5,3'-triiodothyronine). Essential for providing the brain with appropriate levels of T3 during the critical period of development.

Gene ID:

DIO2

Uniprot

Q92813

Synonyms:

Type II iodothyronine deiodinase (EC 1.21.99.4) (5DII) (DIOII) (Type 2 DI) (Type-II 5'-deiodinase), DIO2, ITDI2 TXDI2

Immunogen:

Recombinant Human Type II iodothyronine deiodinase protein (74-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A