

# MTBP Antibody, Biotin conjugated



PACO58439

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Inhibits cell migration in vitro and suppresses the invasive behavior of tumor cells. May play a role in MDM2-dependent p53/TP53 homeostasis in unstressed cells. Inhibits autoubiquitination of MDM2, thereby enhancing MDM2 stability. This promotes MDM2-mediated ubiquitination of p53/TP53 and its subsequent degradation.

**Gene ID:**

MTBP

**Uniprot**

Q96DY7

**Synonyms:**

Mdm2-binding protein (hMTBP), MTBP

**Immunogen:**

Recombinant Human Mdm2-binding protein (567-724AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A