

KIFC1 Antibody, HRP conjugated



PACO58505

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Minus end-directed microtubule-dependent motor required for bipolar spindle formation. May contribute to movement of early endocytic vesicles. Regulates cilium formation and structure.

Gene ID:

KIFC1

Uniprot

Q9BW19

Synonyms:

Kinesin-like protein KIFC1 (Kinesin-like protein 2) (Kinesin-related protein HSET), KIFC1, HSET KNSL2

Immunogen:

Recombinant Human Kinesin-like protein KIFC1 protein (337-498AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A