

PTTG2 Antibody



PACO58540

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

nucleus, peptidase inhibitor activity, homologous chromosome segregation, negative regulation of mitotic sister chromatid separation, regulation of transcription, DNA-templated.

Gene ID:

PTTG2

Uniprot

Q9NZH5

Synonyms:

Securin-2 (Pituitary tumor-transforming gene 2 protein), PTTG2


Immunogen:

Recombinant Human Securin-2 protein (1-202AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescence staining of PC-3 cells with PACO58540 at 1:166, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).